ROHO AIRLITE CUSHION



SKU: AL-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The AirLITE cushion is a non-adjustable, maintenance-free foam cushion with a ROHO air insert that conforms to your body while also preventing you from bottoming out in your cushion. The gentle contour of the leg trough helps you maintain the right position and posture, and lets you feel more secure.

FEATURES

- Designed for clients with tactile sensation who are at low risk of developing skin breakdown
- Contoured foam base with an internal sealed air cell provides comfort and support for basic positioning and pressure management needs.
- **“Bottom-Out” Protection:** The non-adjustable ROHO AIR FLOATATION™ air insert built into the cushion under the ischial has been enhanced to provide better bottom-out protection.
- **Softer Foam:** Softer foam allows the cushion to conform to the individual’s physiology and provide a more comfortable seating surface.
- **Durability and Simplicity:** The 24 month warranty shows our confidence in the quality and durability of the product. The AirLITE is an off the shelf, ready to use product. No adjustment is required.
- **Heavy Duty Cover:** The cover design, with full zipper coverage, provides enhanced fluid resistance to protect the foam. At the same time, the cover is 100% PVC free, reinforcing ROHO’s commitment to more environmentally friendly products.
- All models include: 2 way stretch cover, Hand pump & Repair kit, Operating instructions.

SPECIFICATIONS

Product	Total Height	Safe Working Load	Pressure Care Rating	Material
ROHO AirLite	2 - 3.75 in	113 kg	Low Risk	Air Cells

View the ROHO Cushions Selection Guide

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au