

ROYALE MEDICAL HEAVY DUTY BATH TRANSFER BENCH

SKU: S11190

PRODUCT DESCRIPTION

Purpose & Benefits of Royale Medical Heavy Duty Bath Transfer Bench

The Royale Medical Heavy Duty Bath Transfer Bench is engineered to provide a safe, comfortable, and stable platform for individuals with mobility issues to transition into and out of the bathtub. This bench is especially beneficial for bariatric patients and those who find it difficult to stand in the shower, enhancing their ability to maintain hygiene independently in a secure manner.

Features

- **Reinforced powder-coated frame** – supports up to 280 kg for reliable, long-term use
- **Extra-wide seat** – provides comfort and stability for larger body shapes
- **Height-adjustable legs** – accommodate individual needs and ensure a level fit
- **Reversible backrest** – supports left- or right-hand bath access
- **Drainage holes** – reduce water pooling and promote hygiene
- **Suction-tipped leg** – enhances stability inside the bath
- **Tool-free assembly** – quick to install and reposition

Specifications

- **Safe Working Load**: 280 kg
- **Dimensions**:
  - Overall Height: 82.5 cm - 97.5 cm
  - Overall Length: 52 cm
  - Overall Width: 78.5 cm
  - Seat Height: 49 cm - 64 cm
  - Seat Depth: 41 cm

- Seat Width: 68 cm
- Armrest Height: 27.5 cm
- Unit Weight: 6 kg

The Royale Medical Heavy Duty Bath Transfer Bench is an essential aid for individuals who require support while bathing. It provides a stable and secure environment, reducing the risk of slips and falls in the bathroom, thereby promoting independence and preserving dignity.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au