

SHEAR PRESSURE CARE ASSISTANT XD1900

SKU: XDPCA

PRODUCT DESCRIPTION

Purpose & Benefits

The **Shear Comfort Pressure Care Assistant** is designed to provide targeted relief and comfort for specific body areas—such as heels, hips, knees, shoulders, or head—while lying in bed. It securely attaches via linen flaps tucked under the mattress, keeping the wool protection exactly where it's needed. Made from high-density XD1900 medical-grade wool, it offers superior pressure redistribution, a healthy skin microclimate, and long-lasting durability.

Features

- Linen flaps enable easy secure attachment under the mattress for precise placement
- Constructed with high-density XD1900 wool (1,900 g/m²) for effective pressure redistribution
- Promotes a healthy microclimate—wool absorbs up to 30 % of its weight in moisture vapor without feeling damp, allowing heat and moisture to dissipate
- Machine washable up to 80 °C thanks to fabric-backed construction for hygienic cleaning
- Guaranteed durability for up to 50 washes when cared for according to manufacturer’s instructions

Specifications

- **Brand:** Shear Comfort (Healthcare Innovations Australia)
- **Product Type:** Pressure Care Assistant (bed overlay)
- **Material:** XD1900 medical-grade wool with fabric backing
- **Colour:** Blue
- **Dimensions:** 70 cm × 30 cm
- **Washability:** Machine washable up to 80 °C; up to 50 washes guaranteed

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au