

SHOWERBUDDY PETITEBUDDY

SKU: P1-H

PRODUCT DESCRIPTION

Purpose & Benefits of the ShowerBuddy PetiteBuddy

The ShowerBuddy PetiteBuddy is an accessory designed to adapt standard ShowerBuddy transfer systems for smaller users, including children and petite adults. By modifying the chair's dimensions, the PetiteBuddy ensures a comfortable and secure fit, enhancing safety and usability during bathroom transfers. This adaptability allows the ShowerBuddy system to grow with the user, providing a cost-effective solution that eliminates the need for frequent equipment upgrades.

Features

- **Adjustable Backrest and Neck Pillow**: Offers customisable support to accommodate the user's specific comfort and positioning needs.
- **Reduced Width Between Armrests**: Narrows the space between armrests to better fit smaller body frames, ensuring a secure and comfortable seating position.
- **Smaller Commode Opening**: Features a reduced commode opening size, providing added safety and comfort for petite users during transfers.
- **Tool-Free Assembly**: Designed for easy installation without the need for tools, allowing for quick setup and adjustments as needed.
- **Compatibility**: The PetiteBuddy accessory is compatible with ShowerBuddy models SB1, SB2, SB2T, and SB3T, offering versatility across various transfer systems.

Specifications

- **Backrest Width**: 380 mm
- **Backrest Depth Adjustment**: 325 mm to 375 mm
- **Width Between Armrests**: 420 mm

- **Commode Opening Size**: 170 mm
- **Compatibility**: Compatible with ShowerBuddy models SB1, SB2, SB2T, and SB3T
- **Assembly**: Tool-free assembly for easy installation

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au