


# THERAMED COCCYX TAILBONE WEDGE

**SKU:** A140412001-H

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

The Tailbone Wedge Seat Cushion relieves coccyx pressure.It has a coccyx cut-out designed to distribute weight evenly and relieve painful coccyx pressure.

## The Coccyx Tailbone Wedge Cushion is designed to alleviate a range of conditions:

- Relief for pain or discomfort in the coccyx region
- Lower back pain during pregnancy
- Poor posture accompanied with coccyx pain

## TWO COVER OPTIONS AVAILABLE

- SteriPlus: Moisture and bacteria resistant, features a convenient zip for easy removal, and can be wiped down with a cloth or removed and put in the washing machine.
- DuraFab: Attractive, soft, durable and washable polyester fabric.

## SPECIFICATIONS

| | |
|---|---|
| **Medium Size** | 42cm across x 40 cm x 2.6 cm (thin) x 7 cm (thick) |
| **Inner** | 100% open cell 'breathing' urethane foam. |
| **Dura-Fab Outer Cover** | Woven polyester. It is a super-soft easy-clean washable fabric. |
| **Steri-Plus Outer Cover** | Moisture and bacteria-resistant polyurethane ply bonded to a knitted polyester backing. |

**Country of Manufacture** Australian Made

Buy this Coccyx Wedge Cushion by TheraMed online at Medimart, with AfterPay and the NDIS available. Shop Now!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

**Elizabeth Vale Showroom**
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

**Reynella Showroom**
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

## VIC

**Frankston Showroom**
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

**Monash Showroom**
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

## NT

**Darwin Showroom**
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au