

THRONE BARIATRIC TOILET RAIL

SKU: RB150

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Bariatric Toilet Rail is designed to assist the larger build. This particular rail can hold weight up to 150 kg to 180 kg.

With its freestanding feet and rubber grips, the Bariatric rail relies on the floor of the bathroom/toilet, rather than being supported by the toilet itself, providing excellent stability. The legs and feet are set back to be of no hindrance to the user.

SPECIFICATIONS

- Distance Between Handrails: Approx 600 mm
- Height Above Seat: 250 mm (with adjustable legs)
- Overall Width: 750 mm

https://youtu.be/ti6wiFWmTwQ

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au