VARILITE 2 POINT HIP BELT



SKU: 8321-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

VARILITE 2-Point Hip Belt attaches to the wheelchair at two points, with padding placed relative to the ASIS (anterior superior iliac spine) to provide pelvic support or stability. Attach the 2-point hip belt to the side wheelchair canes to promote anterior pelvic tilt, and attach to the back canes, or to a solid back, several inches above the seat-to-back junction for posterior pelvic tilt.

This dual attachment system not only offers customization to suit the user's specific physical needs and seating position but also contributes to increased overall safety by preventing slippage and maintaining optimal pelvic alignment.

FEATURES

- 40 mm webbing with options for different buckles and padding sizes (medium and large)
- The design allows the padding to be positioned to fit each individual and secured in place with tri-glides (attaching hardware)
- Contoured shape laminated foam and binded edges
- To accommodate growth, replace padding with a larger size
- Can be modified to a 4-point hip belt
- Available sizes: M, L

SPECIFICATIONS

Description	Buckle	Closure	Webbing	Pad Dimension	Length
2-Point, Padded Hip Belt, Medium	Push Button, Single Pull	Metal-to-Metal	1 1/2 in. (3.8 cm)	9 × 2 1/2 in. (22 × 6.3 cm)	60 in. (152 cm)

Description	Buckle	Closure	Webbing	Pad Dimension	Length
2-Point, Padded Hip Belt, Large	Push Button, Single Pull	Metal-to-Metal	1 1/2 in. (3.8 cm)	10 5/8 × 2 1/2 in. (27 × 6.3 cm)	60 in. (152 cm)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au