

REDGUM COMPACT SIDE FOLDING WALKER - RED OR BLUE

SKU: RG4401

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This wheeled walker provides the user a safe walking aid and promotes confidence whilst in its operation. The walker comes with a seat for the individual to use when feeling tied. Comes with a vinyl carry bag.

FEATURES

- Ergonomically designed
- Single hand, side folding action for easy storage
- Lightweight aluminium construction and powder coated finish
- Soft operation hand brake
- Seat height at 55cm from floor
- 11.5cm wide padded backrest
- Removable vinyl bag fitted
- Walking stick holder
- Walking stick holder
- 20cm castor wheels

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au